



Presentation:
0101026001 - 12 UN / one (1) 2g syringe

DIAMOND UNIVERSAL

It is a universal polishing paste for use exclusively by dentists or prosthetics.

- INDICATED FOR POLISHING:**
Resin composites, porcelain, dental enamel and other restorative materials.
- EXTRA FINE GRANULATION**
Produced with micronized diamond particles.
- GREAT AFFINITY WITH WATER**
Easy removal after use.
- EXCELLENT PERFORMANCE**
High brightness capacity.



Presentation:
010201003 - 20 UN / 1 25g jar
010201004 - 20 UN / 1 25g jar (with fluoride)

OBTUR

Temporary restorer used to seal the pulp cavity, preventing microleakage and the exit of intracanal medications, ensuring the success of endodontic treatment.

- EXCELLENT MARGINAL SEALING.**
- GOOD CONSISTENCY.**
- OBTUR WITH FLUORIDE HELPS REDUCE THE SENSITIVITY OF THE TREATED SITE**
- FLUORIDE VERSION**

MORIZ TECHNOLOGY

STM MEDITECH
TC/26/1245-3, 2nd Floor, Monarch Building
Palakkal Angadi, Thrissur
Kerala India - 680 001.
E-mail: stmmeditech@gmail.com
Ph: 9946574883, 8590070532

www.morizdental.com

92 | TEMPORARY RESTORER



Presentations:
0102021001 - 20 UN / Blue Abracec
0102021002 - 20 UN / Pink Abracec

ABRITEC

It is a device to assist in opening the patient's mouth, facilitating the execution of clinical procedures.

- SUPER RESISTANT:**
Non-toxic silicone and pigment.
- AUTOCLAVABLE, 121° C.**
- AVAILABLE IN:**
PINK and BLUE.
- GREATER SAFETY:**
Presence of a hole for flossing.

CE



Presentations:
010201003 - 10 UN / One (1) adult: retractor U
010201004 - 10 UN / One (1) children retractor U
010201007 - 10 UN / One (1) adult: retractor V
010201008 - 10 UN / One (1) children retractor V

RETRACTOR U V

Lateral retractor for lips and cheek.

- Retractor in U**
Indicated to aid in intraoral photographs of the front and anterior teeth.
- Retractor in V**
Indicated to aid in lateral intraoral photographs and of the posterior teeth.
- HIGH PERFORMANCE**
It allows great amplitude for frontal photographing from 2nd molar to 2nd molar.
- AUTOCLAVABLE 121° C.**
- AVAILABLE IN THE SIZES:**
ADULT and CHILDREN.

CE

PARAFFIN-BASED WAX FOR BRACKET PROTECTION

AESTHETICS (transparent material) **MENTOLATED SCENT,** non-toxic.

PRESENTATION:
0105001001 - 50 UN
Ten (10) cases with five (5) sticks each

PLASTIC DEVICE TO ASSIST IN THE PASSAGE OF DENTAL FLOSS UNDER ORTHODONTIC FRAMES AND FIXED BRIDGES.

PRESENTATIONS
010201001 - 50 UN.
Package with fifty (50) conductors
010201002 - 50 UN.
Package with ten (10) conductors

THERMOPLASTIC RESIN BOX FOR STORING INTRA AND EXTRAORAL ELASTIC BANDS.

AVAILABLE IN:

PRESENTATION:
0102017003 - SOUND Ten (10) units per package

RETAINER BOX
Box for storing the retainer.

PRODUCED IN THERMOPLASTIC RESIN.
great resistance.

AROMATIZED BOXES TUTTI-FRUTTI.

Pink Orange Green Yellow Blue
Pearly blue Silver Pearly white

PRESENTATION:
0102007008 - 50 UN / Multicolor
0102007002 - 16 UN / Pearly blue
0102007005 - 16 UN / Pearly white
0102007007 - 16 UN / Silver
0102007009 - 16 UN / Pink

PACKAGING FOR ORTHODONTIC ELASTIC BANDS

Packaging for storage and organization of orthodontic elastic bands.

PRESENTATION:
0104018001 - 20 UN
Fifty (50) package

www.maquira.com.br

FACIAL MASK FOR PROTRACTION OF THE JAW (REVERSE TRACTION).

PRESENTATION:
One (1) facial mask + two (2) adhesives.
0102023001 - 12 UN / Blue
0102023002 - 12 UN / Pink
0102023002 - 12 UN / Gray



Presentations:
0101025001 - 10 UN / one (1) 2g syringe (white color)
0101025002 - 10 UN / one (1) 2g syringe (shaded color)

Presentation:
0101025003 - 10 UN Five (5) 2g syringes (opaque white) + 1.37%
Acid gel
0101025004 - 10 UN Five (5) 2g syringes (shaded color) + 1.37%
Acid gel

MAX SEAL

Light-curing resin sealant for use in pits and fissures of posterior and permanent teeth with complex anatomy to prevent caries.

- INDICATED FOR PATIENTS WITH HIGH CARIOGENIC ACTIVITY AND ON NEWLY ERUPTED TEETH.**
- MAKES A RESISTANT BARRIER**, which hinders the adhesion of plaque and food remains.
- CONTAINS TWO SOURCES OF FLUORIDE:** sodium fluoride and calcium fluoride.
- EXCELLENT CONSISTENCY**, it does not run down from the site.
- PHOTOPOLYMERIZABLE:** 20 seconds.
- READY FOR USE**, application is made directly from the syringe, with no need for manipulation.
- AVAILABLE IN:** SHADED OPAQUE WHITE.



Presentations:
0107003001 - 10 UN / #500 / One (1) 250cm Rask
0107003002 - 10 UN / #600 / One (1) 250cm Rask
0107003003 - 10 UN / #0 / One (1) 250cm Rask
0107003004 - 10 UN / #1 / One (1) 250cm Rask
0107003005 - 10 UN / #2 / One (1) 250cm Rask

RETRACTOR

Cord for gingival retraction

- CORD THICKNESSES:**
000 Ultra Extra Thin
00 Extra Thin
0 Thin
1 Medium
2 Reg. Thick
- 100% COTTON CORDS**
- KNITTED CORDS:**
High absorptive power of hemostatic materials and also so that the preparation remains free of moisture.
- FINE CORD**, allows easy positioning and remain fixed as desired.
- PACKING WITH CUTTING AREA**, more practicality.



Presentations:
0101022001 - 12 UN / one (1) 2g syringe blue
0101022003 - 12 UN / one (1) 2g syringe green

MAXDAM

Indicated for protection of gingival tissue in cases of risk of irritation by contact with products used in the clinical procedure, such as, for example, tooth whitening with Hydrogen Peroxide or Carbamide.

- DOES NOT IRRITATE GINGIVAL TISSUES**
- EXCELLENT SEALING**, does not move from the applied location.
- LIGHT-CURING RESIN:** 30 seconds.
- EASY TO REMOVE:** it comes out entirely at the time of removal.
- AVAILABLE IN:** BLUE and GREEN.



Presentations:
0104008002 - 5 UN / KIT with 50 matrices, 8 protectors and 2 clips

Presentations:
0104008003 - 5 UN / KIT with 25 matrices, 4 protectors and 1 clip

UNIVERSAL MATRIX

System composed of preformed and disposable sectional matrices and clips for fixation.

- INDICATED FOR CLASS II RESTORATIONS.**
- BEST ACCOMMODATION IN THE CAVITY.**
The design of matrices favors the restoration of dental anatomy.



Presentations / 5 UN:
0104008004 - 50un refill S
0104008005 - 50un refill M
0104008006 - 50un refill L
0104008007 - 50un misc.

Presentations / 5 UN:
0104008008 - 2 clips + 8 protectors
0104008009 - 20 protectors





Presentations:
0102019001 - 30 UN / one (1) blue dappen dish
0102019002 - 30 UN / one (1) green dappen dish
0102019004 - 30 UN / one (1) pink dappen dish

SILICONE DAPPEN DISH

Suitable for handling and storing products for dental, medical and pharmaceutical applications.

INDICATED FOR PRODUCTS:
Composite, prophylactic
paste and others.

3 CAVITIES
for handling.

AVAILABLE IN
BLUE, PINK and GREEN.

AUTOCLAVABLE,
121° C.



Presentations:
0102015001 - 20 UN / One (1) green dish
0102015002 - 20 UN / One (1) blue dish

AMALGAM HOLDER

Indicated for the transport and insertion of amalgam in the dental cavity.

MATERIAL
PLASTIC

AVAILABLE IN
PINK and GREEN.

AUTOCLAVABLE,
121° C.

CE



PLASTIC IMPRESSION TRAY

They are used to obtain impression molds for teeth, using alginates or any other impression material.

PRESENTS
ROUNDED AND
SMOOTH EDGES.

VERTICAL
DRILLING, increases
material retention.

AUTOCLAVABLE
121° C

CE

Presentation:
0102005001 - 10 UN / Tray No. 1
0102005002 - 10 UN / Tray No. 2
0102005003 - 10 UN / Tray No. 3
0102005004 - 10 UN / Tray No. 4
0102005005 - 10 UN / Tray No. 5
0102005006 - 10 UN / Tray No. 6
0102005007 - 10 UN / Tray No. 7
0102005008 - 10 UN / Tray No. 8
0102005009 - 4 UN / Tray KIT 1 to 8



Presentations:
0102010001 - 20 UN / Twelve (12) strips 125 x 4.0 x 0.12 mm
0102010002 - 20 UN / Twelve (12) strips 125 x 6.0 x 0.12 mm

AIRON STEEL ABRASIVE STRIP

Steel strip for finishing and polishing the interproximal faces of dental restorations.

INDICATED FOR:
Composite resin, acrylic
resin, glass ionomer.

EXCELLENT
INTERPROXIMAL WEAR
(between the teeth).

SIZES:
125 x 4 mm
125 x 6 mm

CE



Presentations:
0102010001 - 50 UN / Fifty (50) strips.

AIRON POLYESTER STRIP

Indicated to assist in dental restorations with composite resin and to isolate adjacent teeth.

100% POLYESTER, SUPER
RESISTANT.

TRANSPARENT STRIP,
improves the passage of
light during light curing of the
restorative material.

SIZES:
120 x 10 mm

CE



Presentations:
0102010001 - 20 UN / One hundred and fifty (150) strips

AIRON POLYESTER ABRASIVE STRIP

Interproximal finishing and polishing strips for composite resin dental restorations.

FLEXIBLE STRIPS,
adapt to the contour of
proximal regions.

NEUTRAL CENTER,
Facilitates insertion between
teeth.

2 WEAR THICKNESS,
a thin part and a thick part,
allows wearing out the region
and also finishing and polishing

SIZES:
170 x 3.8 mm

CE



Presentations:
0102001005 - 10 UN / Chemical adult
0102001006 - 10 UN / Chemical children
0102001011 - 10 UN / Autoclavable adult
0102001012 - 10 UN / Autoclavable children

EXPANDEX RETRACTOR

Lateral retractor for lips and cheek.

FOR PROCEDURES:
Bracket bonding
Intraoral photography
Application of fluoride
Adhesive restorations
Prophylaxis
Whitening
Desensitizers
Prosthetic preparations.

AVAILABLE IN VERSIONS:
Autoclavable 121° C.
Chemical Sterilization

AVAILABLE IN THE SIZES:
ADULT and CHILDREN.

CE



Presentations:
0102001014 - 10 UN / Large retractor
0102001015 - 10 UN / Medium retractor

FLEX RETRACTOR

Lateral retractor for lips and cheek.

INDICATED FOR:
Assistance in procedures that require
retreat of lips and cheeks for greater
dental exposure and without the
interference of the patient's tongue,
among them: dental whitening in the
office, application of desensitizers,
bonding brackets, prophylaxis, application
of adhesive restorations, minor
anterior gingival surgery and prosthetic
preparations.

WITH BITE REST:
Greater convenience for the patient.

AUTOCLAVABLE 121° C.

AVAILABLE IN THE SIZES:
ADULT and CHILDREN.

CE



Presentation:
0102000307 - 12 UN / Twenty (10) aspirators

ENDO FLEX ASPIRATOR

Endodontic suction device used to assist in the intracanal suction of blood, fluids and irrigation solutions.

- MULTIFUNCTIONAL TIP,** can be adjusted to the desired angle.
- EXTREMELY FINE AND FLEXIBLE TIP,** easy suction in difficult-to-access root canals.
- PACKAGE STERILIZED INDIVIDUALLY** (ethylene oxide).
- PERFECT FIT ON THE HOSE.**

73 | ENDODONTICS



Presentation:
0104004001 - 50 UN / Twelve (12) super carbon basic sheets
0104004002 - 50 UN / Twelve (12) super carbon sheets

SUPER CARBON

To register occlusal contacts in the procedures for adjusting restorations, prosthetic parts and dental surfaces.

- RESISTANT:** Resistant to traction and moisture.
- BLUE and RED,** Double face super carbon.
- 100µm THICKNESS** not interfering with occlusion.
- SIZES:** Super carbon - 8.7 x 1.8 mm
Basic - 10.3 x 1.8 mm



Presentation:
0102016001 - 30 UN / one (1) blue dappen dish
0102016002 - 30 UN / one (1) pink dappen dish
0102016004 - 30 UN / one (1) green dappen dish

PLASTIC DAPPEN DISH

Suitable for handling and storing products for dental, medical and pharmaceutical applications.

- INDICATED FOR PRODUCTS:** Composite, acrylic, prophylactic paste and others.
- AVAILABLE IN** BLUE, PINK and GREEN;
- AUTOCLAVABLE,** 121° C.

41 | ACCESSORIES



Presentation:
0102013001 - 30 UN / one (1) blue plastic "tamborel"
0102013004 - 30 UN / one (1) pink plastic "tamborel"
0104010001 - 30 UN / one (1) aluminum "tamborel"
0104009001 - 30 UN / "tamborel" refill with 50

TAMBOREL

Plastic or aluminum support to hold endodontic files.

- ACCOMMODATES THE LARGEST NUMBER OF** ENDODONTIC LIMES.
- AVAILABLE IN:** BLUE and PINK.
- "TAMBOREL" REFILL 100% POLYESTER**
- AUTOCLAVABLE 121°C**



Presentation:
0104011001 - 50 UN / one (1) aluminum ruler
0102006002 - 50 UN / one (1) plastic ruler

ENDODONTIC RULER

Millimeter ruler (plastic and aluminum) for measuring endodontic files.

- SUPER RESISTANT RULER,** resistant to multiple sterilization cycles.
- PLASTIC RULER SIZE:** 35mm.
- ALUMINUM RULER SIZE:** 40mm.

CE

72 | ENDODONTICS



Presentation:
0101034001 - 20 UN / One (1) 15g powder flask + one (1) 10g liquid flask

IONGLASS C

Conventional self-curing glass ionomer cement for cementation.

- INDICATED FOR:** Fixation of Crowns
Fixation of Prostheses
Temporary cementation of prosthetic parts
Fixation of orthodontic bands.
- GREATER CAPACITY TO RELEASE FLUORIDE.**
- RADIOPAQUE.**
- HIGH CAPACITY OF ADHESION TO DENTAL STRUCTURE.**



Presentation:
101034004 - 20 UN / A2
0101034007 - 20 UN / A3
One (1) 15g powder flask + one (1) 10g liquid flask

IONGLASS R

Conventional self-curing glass ionomer cement for restoration.

- INDICATED FOR:** Restoration of deciduous teeth, Class III and V restorations, repair of erosions in cervical regions. For atraumatic restorative treatment (ART). Provisional restorations.
- QUICK SETTING CEMENT.**
- GREATER CAPACITY TO RELEASE AND REFILLING FLUORIDE.**
- AVAILABLE IN THE COLORS A2 / A3**

58 | SURGICAL



Presentations:
0101033001 - 20 UN / one (1) 4mL. Rack
0101033002 - 20 UN / one (1) 5mL. Rack
0101033003 - 20 UN / one (1) 5mL. Rack

BOND 2.1

Conventional two-step light-curing adhesive system for enamel and dentin.

PRIMER & BOND, in the same bottle.

STRONG ADHESIVE RESISTANCE

ETHANOL-BASED SOLVENT

PHOTOPOLYMERIZATION: 20 seconds.



Presentation:
01010330001 - 12 UN / one (1) 5mL. Rack

SILANO

Chemical bonding agent for adhesive cementation work.

INDICATED FOR: Ceramics, ceromers, composite resins and fiberglass posts.

INCREASE THE ADHESIVE CAPACITY BY UP TO 30%.



Presentations:
0101044001 - 6 UN / One (1) 50g base tube + one (1) 15g accelerator tube + 1 mixing pad

PROVISORY

Cement based on zinc oxide and organic acids for temporary fixation of prosthetic work.

CEMENT BASED ON ZINC OXIDE AND ORGANIC ACIDS

RADIOPAQUE

HIGH RESISTANCE to withstand chewing loads.



Presentations:
0101048002 - 6 UN / one (1) 50g base tube + one (1) 15g accelerator tube + 1 mixing pad

PROVISORY NE

Cement based on zinc oxide and organic acids, WITHOUT EUGENOL for temporary fixation of prosthetic work.

CEMENT BASED ON ZINC OXIDE ORGANIC ACIDS

NO EUGENOL does not interfere with other components

RADIOPAQUE.

HIGH RESISTANCE, to withstand chewing loads.

45 | ADHESIVES AND SILANE



ZIRCONFILL

Nanohybrid composite resin

INDICATED FOR RESTORATION PROCEDURES

ENAMEL AND DENTIN
All classes (I, II, III, IV, and V). Aesthetic corrections, fragment sticking.

DIRECT VENEERS
Made of composite resins. Indirect restorations (inlays and onlays).

CONFECTION OF THE FILLING CORE

THE ONLY COMPOSITE RESIN CONTAINING DIATOMITE

- EXPRESSIVE REDUCTION of polymerization contraction;
- HIGH LEVEL OF CONVERSION from monomer to polymer;
- HIGH superficial smoothness, LESS AGGREGATION of dental biofilm;
- CLINICAL LONGEVITY of color stability.
- HIGH TIXOTROPY and workability.

PRESENTATIONS One (1) 4 gram syringe

0101047001 Sun - DA1	0101047017 Sun - EB1
0101047002 Sun - DA2	0101047018 2un - EB2
0101047003 2un - DA3.5	0101047019 2un - EB3
0101047004 2un - DA3	0101047020 2un - EC2
0101047005 2un - DA4	0101047021 2un - EC3
0101047006 Sun - DB1	0101047022 2un - D-BLEACH H
0101047007 Sun - DB2	0101047023 2un - D-BLEACH L
0101047008 2un - DB3	0101047024 2un - D-BLEACH M
0101047009 2un - DC2	0101047025 Sun - DA2
0101047010 2un - DC3	0101047026 Sun - DA3
0101047011 2un - D-BLEACH	0101047027 2un - T-BLUE
0101047012 Sun - EA1	0101047028 2un - T-GRAY
0101047013 Sun - EA2	0101047029 2un - T-NEUTRAL
0101047014 Sun - EA3.5	0101047030 2un - T-ORANGE
0101047015 Sun - EA3	0101047031 2un - T-YELLOW
0101047016 2un - EA4	

53 | SURGICAL

ESSENTIAL FOR PROTECTION OF PERIODONTAL SURGERY

PERICEM

PERIODONTAL SURGICAL CEMENT

INDICATED FOR:

Providing safe local protection for surgical wounds, keeping them clean. Cement is produced by mixing its components (base paste and accelerator paste) at the time of use.



EXCELLENT PERFORMANCE

PERFECT PROTECTION OF SURGICAL WOUNDS, keeping them clean.

NO EUGENOL, ensuring its use does not irritate the oral mucosa.

FLAVOR PEPPERMINT.

GREAT SETTING TIME: 2 to 3 minutes.

EXCELLENT TEXTURE and plasticity after setting.



PRESENTATION:

ONE (1) 90g BASE PASTETUBE
ONE (1) 90g ACCELERATOR PASTE TUBE
ONE (1) MIXING PAD
0101043001 - 6 UN

HEMOSPON[®] TAPE
Hemostatic sponge

PAGE 60

THERMAL BAG Hot and cold

PAGE 63

SMALL
SIZE, GREAT
POSSIBILITIES

FIBER POST

FIBERGLASS POST

Pin manufactured with a fiberglass- and epoxy resin-based resin matrix for reconstruction of teeth that have undergone endodontic treatment and require restoration.

It acts with intraradicular reinforcement of the dental structure, promotes retention for the definitive restorative material or filling core, in the case of indirect restorations.



SEQUENCE OF PREPARATION WITH LARGO

N° 0,5	N° 1
N° 1	N° 1 + 2
N° 2	N° 1 + 2 + 3
N° 3	N° 1 + 2 + 3 + 4

EXCELLENT PERFORMANCE

- DOUBLE CONICITY** excellent adaptation to root canals.
- GREAT RESISTANCE.**
- Optimum TRANSLUCIDITY**, allowing photopolymerization of resin cements.
- EXCELLENT RADIOPACITY**, better control on radiographs.
- SIZES:** No. 0.5 - 1 - 2 - 3

PRESENTATIONS

Five (5) PINS
0108001002- 10 UN / No. 0,5
0108001003- 10 UN / No. 1
0108001004- 10 UN / No. 2
0108001005- 10 UN / No. 3



DUAL FORCE
Dual resin cement

PAGE 51

SILANE
Coupling agent

PAGE 45

maquira
Dental Group



Presentations:
0101004001 - 24 UN / One (1) 10% acid gel syringe with 2,5mL
0101004002 - 24 UN / One (1) 5% acid gel syringe with 2,5mL

5% & 10% HYDROFLUORIC ACID GEL

Acid conditioner for dental porcelain surfaces, helping to promote retention in restorations.

- 5% INDICATED FOR:** Vitreous ceramics based on lithium disilicate.
- 10% INDICATED FOR:** Feldspathic Porcelain, Vitreous ceramics reinforced with leucite.
- GREAT AFFINITY WITH WATER,** easy to remove after application.
- 5% ACID YELLOW COLOR, 10% ACID RED COLOR** colors facilitating visualization in applied location.

CE



Presentations:
0101005001 - 50 UN / Three (3) 37% acid gel syringes with 2,5mL

37% PHOSPHORIC ACID GEL

Acid conditioner for enamel and/or dentin surface, allowing greater and better retention of resin materials.

- GREAT AFFINITY WITH WATER,** easy to remove after application.
- EXCELLENT VISCOSITY,** remains in place applied.
- BLUE COLOR,** facilitates visualization in place applied.
- APPLICATION PERIOD:** Enamel: 30 to 15 seconds, Dentin 15 seconds

CE

46 | CONDITIONING AGENTS

Composite resins with
HIGH QUALITY
STANDARD

Applic Flow

FLUID COMPOSITE RESIN

INDICATED FOR THE PROCEDURES

- Direct restorations on anterior and posterior teeth;
- Minor restorations;
- Class III and V restorations.;
- Sealing of pits and fissures;
- Lining of cavities;
- Repair in composite resin;
- Repair of enamel defects;
- Dental fragment sticking;
- Tunnel type preparation;
- Planification of preparation walls.

EXCELLENT PERFORMANCE

- EXCELLENT FLUIDITY** in difficult-to-access cavities.
- HIGH MECHANICAL RESISTANCE.**
- DOES NOT RUN DOWN** from the applied site.
- READY TO USE**
- Excellent RADIOPACITY**
- HIGH BRIGHTNESS POWER.**
- STABILITY OF COLOR** and durability.
- COLORS:** A1 / A2 / A3 / A3,5 / B2.

PRESENTATIONS

ONE (1) 2g SYRINGE

0101036001 - 5 UN / A1
0101036002 - 5 UN / A2
0101036003 - 5 UN / A3
0101036004 - 5 UN / A3,5
0101036005 - 5 UN / B2

UNIVERSAL MATRIX

PAGE 39



BOND 2.1
Dental adhesive

PAGE 45



Applic

NANOHYBRID LIGHT CURING COMPOSITE RESIN

INDICATED FOR THE PROCEDURES

- Occlusal, proximal and occlusal-proximal restorations of posterior teeth;
- Direct veneers in composite resin; Teeth fragment sticking;
- Reduction or closing of diastemas;
- Structural corrections: imperfect amelogenesis, enamel hypoplasias, non-carious cervical lesions;
- Restoration of deciduous teeth.

EXCELLENT PERFORMANCE

- High mechanical resistance.**
- Simplified color system.**
- Color stability.**
- Use of glass filler of German origin; particle with an average size of 0.7 µm.**

TESTED AND APPROVED



Charts: Physico-chemical and mechanical characterization of Applic Composite Resin. Source: Quality control, Maquira.

PRESENTATIONS

ONE (1) 4g SYRINGE

- 5 UN / EA1
- 5 UN / EA2
- 5 UN / EA3
- 5 UN / EA3,5
- 5 UN / EA4
- 5 UN / EB1
- 5 UN / EB2
- 5 UN / EB3

- 5 UN / EC2
- 5 UN / DA1
- 5 UN / DA2
- 5 UN / DA3
- 5 UN / DA3,5
- 5 UN / DB2
- 5 UN / INCISAL

DIAMOND UNIVERSAL
Polishing paste

PAGE 32



POWER BLOCK
Oxygen block gel

PAGE 22



RELATED PRODUCTS